

NEW SMARTARM NIBP SIMULATOR UPGRADE AVAILABLE

Version 7.19A

New Neonatal Software Adds Widest Range of Simulated Neonatal Pressures

This upgrade provides the following enhancements:

- 1. Neonatal Simulation: 6 Neonatal Presets from 35/15 to 150/95**
- 2. Neonatal Pulse Rate Range of 30 to 240 bpm**
- 3. New Simulation Engine: Provides unmatched repeatability and performance**

Customers interested in upgrading their current SmartArms should have the serial number of their SmartArms available when contacting Clinical Dynamics. Customers can contact Clinical Dynamics in three ways:

- 1. Phone: 800-247-6427 in the U.S.A
or 203-269-0090 for international customers**
- 2. Email: service@clinicaldynamics.com**
- 3. Fax: 203-269-3402**

**Clinical Dynamics Corporation
10 Capital Drive
Wallingford, CT 06492
Phone: 800-247-6427
International Phone: 203-269-0090
Fax: 203-269-3402
Web Site: www.clinicaldynamics.com**